

HUGE

Ear Impression Silicone

PERFIT A/RAPID is a low-viscosity A Silicone impression material used to make impressions of the external auditory canal. It is widely used in the production of hearing aids, earplugs and other customized ear products.

Advantages

- **Fast Setting** - Extra short time in the ear (**just 90 sec.!**)
- **Self-Lubricating** - Easy to remove from the ear canal
- **High Thixotropy** - Easy to handle with excellent stability
- **Low Hardness** - Ensures comfortable and safe experience
- **High Elasticity** - Ensures precise and accurate impressions
- **Scannable** - Suitable for digital workflows

Packaging:

- 50ml*8 Cartridges/Box

Technical Parameters

| Technical Parameter | A/RAPID |
|----------------------|---------|
| Working Time (23°C) | 30s |
| Setting Time (23°C) | 90s |
| Final Hardness (24h) | 40 A |
| Color | Blue |
| Shelf Life | 3 Years |

*Date from HUGE Lab

Clinical Procedure

- Cartridge System



Shandon Hude Dental Material Corporation

Add / No.68 Shanhai Road, Donggang District, Rizhao City, Shandong Province,
276800, P.R. China.

Tel / +86 (633) 2277268 marketing@hugedental.com www.hugedental.com

Facebook Huge Dental

Instagram Huge Dental

Youtube Huge Dental